twilight

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Melt Flow Indexer

Operation Manual

The document was written in the technical documentation department of the company LARYEE Technology Co., Ltd.

All rights at this documentation, in particular the rights of the duplication and spreading as well as the translation are at LARYEE Technology, also for the case of patent right registrations. No part of the documentation may be reproduced in any form without previous written agreement of the company LARYEE Technology or processed, multiplied or spread using electronic system. Mistake and technical subject to change.

LARYEE Technology Co., Ltd[®]

LARYEE Technology Co., Ltd is not responsible for any errors in this documentation. A liability for indirect and direct damage, which develops in connection with the supply or the use of this documentation, is impossible, as far as this is legally permissible.

Contents

I.	Instruction.....	2
II.	Operating Principle.....	2
III.	Technical Parameters.....	3
IV.	Structure.....	4
V.	Installation and Adjustment.....	5
VI.	Operation.....	5
VII.	Cautions.....	11
VIII.	Appendix.....	11

I. Instruction

According to the Standard of ISO1133,ASTM1238,ASTM3364, controlled by the PLC, the series of testers are used to, under the high temperature, measure the melt mass-flow rate (MFR) and the melt volume-flow rate (MVR) of thermoplastic, such as polyethylene, polystyrene, polypropylene, ABS resin, melt formaldehyde resin, polycarbonate, nylon, fluorine plastic, and fragrant sulphone. The machine automatically controls the temperature and displays the temperature and the time on a build-in LCD. Two types of cut methods are available: time controlling cut and setting displacement cut. The tester also can automatically calculate MVR value, draw MVR curves and print out all the required experimental datas.

II. Operating Principle

2.1 The tester melts the determined specimen under high temperature by heating the furnace. The melted specimen flows out a specified hole under the given load. In plastic producing activities, we usually use the melt index to, under melted status, measure the physical properties of the polymer, such as the fluidity, the mucosity and so on.

2.2 The melt mass-flow rate is the weight or the volume of the melt specimen. There are two methods, MFR and MVR.

2.2.1 The mass-flow rate **MFR** is the average mass of substance which passes through a given surface per 10 minutes. The unit is (**g/10min**), Denoted by **MFR**, which is defined as

$$\mathbf{MFR=600\times m\div t}$$

Where

- MFR** The melt mass-flow rate of material, unit **g/10min**;
- m** Average mass, unit **g**;
- t** Time interval, unit **s**;
- 600** Reference time (10 min), unit **s**.

2.2.2 The melt volume-flow rate of specimen is the average volume of substance which passes through a given surface per 10 minutes. The unit is (**cm³/10min**). In addition, we define the function of **MVR** as following

$$\mathbf{MVR=42.7\times L\div t}$$

Where

- MVR** the melt volume-flow rate of material, unit **cm³/10min**;
- L** the distance of the predetermined movement of the piston, unit **mm**;
- t** Time interval, unit **s**;
- 42.7** the reference time (**10 min**) times the area of the piston.

2.2.3 MFR Example

There is a group of plastic samples which have been cut per 30 seconds respectively. The weights are 0.0838 g, 0.0862 g, 0.0815 g, 0.0895 g, 0.0825 g respectively. The average weight, $m=(0.0838+0.0862+0.0815+0.0895+0.0825)\div 5=0.0847$ g. Substituting into the equation, we have

$$\mathbf{MFR=600\times m\div 30=1.694(g/10min)},$$

so that the MFR is 1.694g/10min.

2.2.4 MVR Example

In 135 seconds, the tester extrudes 30mm melted plastic specimen.

Substituting the values into the equation, we have

$$MVR = 42.7 \times \pi \times \frac{1}{4} \times \frac{1}{30} = 42.7 \times 30 \div 135 = 9.489 \text{ (cm}^3/10\text{min)}$$

2.2.5 We can calculate the specimen melt density both via **MFR** and **MVR**. In **MVR** example, weights up the extruded specimen, **m**, unit **g**. The specimen melt density expresses by **D**, so that we have:

$$D = m \div MVR \times 600 \text{ (g/cm}^3\text{)}$$

III. Technical Parameters

3.1 Specification of the specimen extruder

- Inner diameter of the die **Φ2.095±0.005mm**
- Length of the die **8.000±0.025mm**
- Inner diameter of the cylinder **Φ9.550±0.025mm**
- Length of the cylinder **160mm**
- Diameter of the piston head **Φ9.475±0.015mm**
- Length of the piston head **6.350±0.100mm**

3.2 Specification of load

There are 7 Standard test resolution levels loading upon the incremental method.

3.2.1. Level 1

$$0.325\text{kg} = (\text{piston rod} + \text{adiabatic cover} + 1\#\text{weight})\text{kg} = 3.187\text{N}$$

3.2.2. Level 2

$$1.200\text{kg} = (\text{level 1} + 2\#\text{0.875weight})\text{kg} = 11.77\text{N}$$

3.2.3. Level 3

$$2.160\text{kg} = (\text{level 2} + 3\#\text{0.960weight})\text{kg} = 21.18\text{N}$$

3.2.4. Level 4

$$3.800 = (\text{level 3} + 4\#\text{1.640weight})\text{kg} = 37.26\text{N}$$

3.2.5. Level 5

$$5.000\text{kg} = (\text{level 4} + 5\#\text{1.200weight})\text{kg} = 49.03\text{N}$$

3.2.6. Level 6

$$10.000\text{kg} = (\text{level 5} + 6\#\text{5.000weight})\text{kg} = 98.07\text{N}$$

3.2.7 Level 7

$$21.6000\text{kg} = (\text{level 6} + 7\#\text{2.500} + 8\#\text{4.100} + 9\#\text{5.000weight})\text{kg} = 211.82\text{N}$$

3.3 Temperature of the cylinder

3.3.1 Temperature Range Set the temperature randomly in the range of 50°C to 400°C

3.3.2 Reach the specified temperature in a moment by PID.

3.3.3 Resolution of the digital display 0.1°C.

3.3.4 Temperature Fluctuation

Fluctuation at the 10mm position of upper end $\leq \pm 0.5^\circ\text{C}$ (50 to 250°C)

Fluctuation at the 10mm position of upper end $\leq \pm 1^\circ\text{C}$ (250 to 400°C)

3.3.5 Temperature gradient

Temperature gradient at the top 10~150mm position of the extruder is less than 1°C.

3.4 Time display

The time is automatically calculated by the build-in clock, Time tolerance $\pm 1\text{s/h}$

3.5 Power **450W, ~220V±10%, 50Hz**

3.6 Dimension **520mm×270mm×550mm**

IV. Structure

4.1 Control Panel

4.1.1 **Monitor** It is a touch screen, English operating interface are available. Operate the tester according to the relative dialog box. Touch the relative dialog box or button to input the test parameters and the relative page shows the test temperature, the rod displacement and the real time curve of MVR.

4.1.2 **Printer** When the test ends, input the mass value of the melted specimen to check or print the MFR test data as well as check or print the data directly.

4.1.3 **PLC** Receive the orders from the touch screen to measure and control the furnace temperature and determine the displacement of the load rod.

4.2 Furnace Chamber shown in fig.1. It consists of heating cover, cylinder, die, load


Fig 1

rod and removable weights.

4.3 The spatula is made up of the motor, the cover and the edge. The motor drives by the tester and rotates automatically to make the edge cut the sample so as to avoid the human error of manual cut.

Note ! Set a proper distance between the edge and the die to guarantee no conglutinated specimen.

4.4 **Displacement Sensor** The displacement sensor installed at the upper right back of the tester, together with the load rod, measures the displacement of the load rod.

V. Installation and Adjustment

5.1 Installation

5.1.1 The Tester should be installed in the lab that has no dust, no vibration, no strong magnetic field and no strong disturbing signal. Place the tester on the workbench.

5.1.2 Put the bubble into the cylinder and observe it. Adjust the foot to make the tester be horizontal.

Note! The adjustment should be done when the temperature of the cylinder is less than 36 °C.

5.2 Adjustment

5.2.1 Turn the spatula and adjust it to just reach the die, and then turn it to the position that it won't shield the extruder.

5.2.2 Installation and download of the Die

Put the die into the cylinder. Use the bolt to poke it out.

5.2.3 Set the test temperature

Set the test temperature according to the physical property of the determined material. Press 'Heat' to heat and control the temperature automatically, then a dialog box prompts the temperature has been constant.

VI. Operation

6.1 Prepare the material for about 4 grams. Select the test parameters in accordance with the standard or test methods, such as test temperature, load, sampling time or set the sample section number for no more than 8 sections

Note! Generally place the specimen for 24 hours under 23°C before the test. Please refer to the relative standard.

6.2. Place the die in the cylinder.

6.3. Preheat the piston (load rod) in the cylinder.

6.4. Prepare the removable weights for the test.

6.5 Dialog box


6.5.1 The initial interface shows as following, touch the screen to enter the next page.

In this interface,

- | | |
|----|--|
| F1 | Touch it can back to the previous interface. |
| F2 | Adjust the screen light |
| F3 | Set password |
| F4 | Set time |

6.5.2 After touching the screen, we can see

In this interface,	
Temperature range	50~400°C. The maximum test temperature is 400°C
Displacement range	1~30mm, the maximum valid test displacement is 30mm
Load level	7 grade, the maximum load (weight) is 21.6Kg
Show temperature	Show the real time temperature
Time setting	Enter the time setting interface
Displacement setting	Enter the section number setting
Parameters setting	See the 6.5.10
Heat	Start to heat

Note

1. **Time Mode and Displacement Mode are two alternative working modes for the tester. The test result is generated in accordance with the settings.**
2. **The “Heat” button is a multi-functional button. When the heating starts, press the button to stop heating, and vice versa.**
3. **When the temp is up to setting temp,the screen will indicate “Temp,ready”.**


6.5.3 Time Mode interface

6.5.4 Displacement Mode interface

6.5.5 Start to test



6.5.6 Test result interface


Note: On above interface, If pressing the button 'MFR Result' directly, the tester shows MVR test result.

6.5.7 MVR test result interface


Note: If pressing the button of “MFR input the mass of Cut-off”, the tester shows MFR test result.

6.5.8 MFR test result interface


Note: Touch the button ‘Density’ to show the density of each Cut-off.

6.5.9 Density interface


6.5.10 Parameters setting interface (The password is 12345678.You don't need to set the parameters,because we have set it well before deliver it)

Touch the Temp.Calibration and enter the temp.calibration interface as the following picture:

Input the target temperature and use a thermometer measuring the

current temperature, then input the difference between the target temperature and the current in the compensation, then save and finish the setting.

6.6. Set the temperature that you select, and then press 'preheat'. The tester alarms when the temperature is constant. Fill the cylinder after the alarm.

6.7. Material Filling and Loading

After the tester reaches a constant temperature, take the piston down and fill about 4 grams specimen into the cylinder with the filling funnel. Compress it with the packing rod. The filling should be finished in a short time. Put the piston back in the cylinder. After 3.5 to 4 min, load the removable weights on the plate of the load rod.

6.8 Whatever the cut method is, weight the sections up respectively and input the values to obtain the test result and print the result report.

6.9 Take the die out through the packing rod after the test and clean the die, the piston and the cylinder with the cotton cloth. Use the drill to poke out the left material and then clean the hole of the die.

VII. Cautions

7.1 Single phase power must be reliably grounded.

7.2 After confirming the results, the instrument is still at constant temperature; the tester can do the experiment continuously as well as re-input the parameters.

7.3 You should switch off the power if you don't need to use the instrument already.

7.4 The user can't uninstall the instrument optionally.

7.5 The instrument should be in the environment that has no strong electromagnetic disturb.

7.6 The cylinder, the piston, the die should keep clean, no collision, no scoring, the cylinder can't be cleaned by the non-designated tools. After a period of usage, use high-level detergent to clean the outside of cylinder and the piston.

7.7 If you can't take the die out, you can take the cylinder out of the instrument, and then use the bolt to poke the die out.

7.8 The cleaning work is easy to be done under high temperature. But at the moment please be more careful not to get burnt.

7.9 If the cylinder, the piston rod or the die has too much abrasion, the test result will be inaccurate. And you can order the parts in time from our company.

7.10 If the instrument is in disorder, restart the power.

VIII. Appendix

8.1 Dry the specimens with kinds of shapes, such as grain, bar, slice and mould and so forth, before the test in accordance with the plastic type requirements. Weight the specimen referring to the Fig. a as following.

Fig. a Weight of the specimen and time interval of incising

Velocity	Weight	Time interval
(g/10min)	(g)	(s)
0.1 ~ 0.5	3 ~ 4	120 ~ 240

0.5 ~ 1.0	3 ~ 4	60 ~ 120
1.0 ~ 3.5	4 ~ 5	30 ~ 60
3.5 ~ 10	6 ~ 8	10 ~ 30
10 ~ 25	6 ~ 8	5 ~ 10

8.2 Recommend test condition

Choose proper load, temperature and die for the materials below according to Fig. b.

Polyethylene 1, 2, 3, 4, 6

POM 3

Polystyrene 4, 7, 11, 13

ABS 7, 9

Polypropylene 12, 14

Polycarbonate 16

Polyamide 10, 15

Acrylate 8, 11, 13

Cellulose acetate 2, 3

Fig. b Test Conditions

No.	Inner diameter of standard mold, unit: mm	Temperature unit: °C	Coefficient of mold unit: g·mm ²	Load unit: kg
1	1.180	190	146.6	2.160
2	2.095	190	70	0.325
3	2.095	190	464	2.160
4	2.095	190	1073	5.000
5	2.095	190	2146	10.000
6	2.095	190	4635	21.600
7	2.095	200	1073	5.000
8	2.095	200	2146	10.000
9	2.095	220	2146	10.000
10	2.095	230	70	0.325
11	2.095	230	253	1.200
12	2.095	230	464	2.160
13	2.095	230	815	3.800
14	2.095	230	1073	5.000
15	2.095	275	70	0.325
16	2.095	300	253	1.200

8.3 Result of weighting

Cool the specimen after cutting. Place the specimens on the balance (Require an Analytical Balance with the 0.001g precision) and weight them respectively.

Calculate the result according to the formula below:

$$\text{MFR} = 600 \times W/t \quad (\text{g}/10\text{min})$$

Where:

MFR: Melt flow rate, unit: g/min

W: Average weight of specimens

t: Time interval of incising, unit: s

Take two decimal points after getting the test result.

Note! Each specimen should be tested horizontally for two times and calculate the MFR value. If the offset is beyond 10% between the two MFR values, find out the reason.

8.4 Calculation of Melt Flow Density

8.4.1 Choose MVR method under automatic mode and press the button 'parameter' and lighten the indicator of volume, input the value of L: 1~30mm (Optional according to the melt index. Generally, if $\text{MFR} < 1$, input 5mm; if $1 < \text{MFR} < 5$, input 10mm; if $\text{MFR} > 5$, input 30mm).

8.4.2 Weight the specimen in 9.4.1 and record the weight, m. After inputting the parameters, place the specimen and test it. Finally, we have the MVR. So that Melt Flow Density is defined as:

$$D_r = \frac{m_L}{0.7163 \times L} \times 10^3 \quad (\text{g}/\text{cm}^3)$$

Melt Flow Indexer Packing List

Model: XNR-400D

Serial No.:

Category	No	Name	Parameter	Quantity	Note
Product	1	Host machine	set	1	
Technical files	1	Operation manual	copy	1	
	2	Certification	copy	1	
	3	Packing list	copy	1	
Accessories	1	Die	pc	1	in the cylinder
	2	Power line	pc	1	
	3	Scrapings rod	pc	1	
	4	Load rod	pc	1	Also as the Weight of 325
	5	spirit level	pc	1	
	6	Ingate drill	pc	1	
	7	Take-up stock pan	pc	1	
	8	Funnel	pc	1	
	9	Forceps	pc	1	
	10	Anchor bolts	pc	4	

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