

twilight

INSTRUMENTOS DE MEDICIÓN INDUSTRIAL

Base para Viscosímetros Ford Tipo Tripie, Biuged BL-BGD131

www.twilight.mx

 / [twilightsadecv](https://www.facebook.com/twilightsadecv)

 / [twilightsadecv](https://twitter.com/twilightsadecv)

 / [twilightsadecv](https://www.youtube.com/twilightsadecv)

Stand For Flow Cup

It is important for any flow cup to be level for the flow cup with its sample during testing of sample viscosity. Biuged instruments offers a simple type stand. It is made up of three adjustable aluminum pole, and glass plate with a spirit level .

Ordering information:

[BGD 130 -----Square Stand](#)

[BGD 131-----Tripod Stand](#)

[BGD 132---- Stainless Steel Stand](#)

[BGD 133---- Tripod Stand with double-walled thermo jacket](#)



BGD 130



BGD 131



BGD 132



BGD 133

twilight

INSTRUMENTOS DE MEDICIÓN INDUSTRIAL



LLÁMANOS

+52(81) 8115-1400 / +52 (81) 8173-4300

LADA Sin Costo:
01 800 087 43 75

E-mail:
ventas@twilight.mx

www.twilight.mx



/ [twightsadecv](https://www.facebook.com/twightsadecv)



/ [twightsadecv](https://twitter.com/twightsadecv)



/ [twightsadecv](https://www.youtube.com/twightsadecv)